

RADIOSYNOVIORTHESIS:

A new therapeutic and diagnostic tool for canine joint inflammation

John M. Donecker, VMD, MS

Nigel R. Stevenson, PhD

Exubrion Therapeutics, Inc.

Buford, GA 30518

DEFINITIONS

Note: First usages of the following definitions in the text are indicated by **bold type**.

beta particle: A high-energy electron emitted from the nucleus of a radioactive atom; beta particles typically have a wide tissue penetration range of 50–5,000 μm that diminishes over distance, making uniform dosing difficult and possibly necessitating shielding during transport and handling.

colloid: A mixture of insoluble microparticles (particles 1 to 1,000 μm) that remain distributed in solution without precipitating or settling to the bottom; nontoxic colloids are used for binding radionuclides to prevent them from escaping the intra-articular space into systemic distribution.

conversion electron: A low-energy electron released from an atomic shell as a result of radioactive decay, resulting when gamma radiation emitted by the nucleus is transferred to the electron; conversion electrons are monoenergetic in contrast to beta particles.

homogeneous tin-117m colloid (HTC): A novel preparation of the radionuclide tin-117m suspended in a colloid, HTC is well suited for intra-articular administration to treat synovial inflammation caused by traumatic injury, OA, and other arthritides.

radiocolloid: A radionuclide-labeled colloid suitable for intra-articular injection.

radionuclide: An unstable isotope of an atom that emits radiation released from the atomic nucleus. Some radionuclides exist naturally, but those with research and therapeutic applications are usually produced artificially. A radioisotope.

radiosynovectomy (RSV): Refers to removal of the synovium and its replacement with fibrotic tissue, for example when

joints are injected with beta-emitting radionuclides with a wide tissue penetration range. RSV is sometimes used as a synonym for radiosynoviorthesis (RSO).

radio-synovi-orthesis (RSO): Injection into the synovial space of a radioisotope to treat joint inflammation and mitigate chondromalacia when systemic or other traditional therapies have failed to produce a satisfactory response. The goal of RSO is reduction of pain and synovial hypertrophy.

synovectomy: Destruction or surgical removal of the membrane (synovium) that lines an articular joint. Open surgical, chemical, radiation, fulguration, and arthroscopic synovectomies are all options for removing potentially damaging synovium from articular joints.

synoviorthesis: A medical therapy using intra-articular injection of a compound that diminishes the degree of synovial hypertrophy, thereby mitigating pain and the development of inflammation and arthritis. Can be performed by chemical synoviorthesis or radiosynoviorthesis, with the latter being preferred when a suitable radionuclide is available.

tin-117m (Sn-117m): An artificially produced radionuclide of tin with medical applications for localized treatment and imaging. Tin-117m has a half-life of 14 days. Two principal forms of the energy that it emits are (1) conversion electrons that have a short penetration range in tissue (~300 μm), and (2) imageable gamma radiation, which enables monitoring of local distribution in tissue. Tin-117m is metastable, indicated by the “m” suffix, meaning that it is a radioisotope with an energetic nucleus and a relatively long half-life and therefore distinct from highly unstable radionuclides with shorter half-lives.

Synovial inflammation is strongly implicated in the pathogenesis of osteoarthritis (OA) and other arthropathies. Synovitis is a common feature of symptomatic but pre-radiographic OA. This indicates that chronic, early-stage joint inflammation occurs well before significant radiographic changes and drives progression toward cartilage loss, osteophyte formation, bone remodeling, and joint space narrowing.¹⁻³ The pathology of early-onset synovitis and its role in OA has been well characterized. Whereas normal synovium is two or three cell layers thick and devoid of inflammatory cells, synovitis results in a number of profound changes in synovial tissue and the joint micro-environment. These include marked hyperplasia and permeability of the synovial lining, significant over-expression of pro-inflammatory mediators and cytokines, infiltration of inflammatory cells, production of degradative enzymes, synovial neovascularization, and increased serum C-reactive protein, a biomarker of inflammation.^{1,2,4-6} The inflammatory response sensitizes peripheral neurons in synovial tissue, resulting in a pain response.⁷

OA should not be thought of as a single entity but as a sequence of events beginning with joint injury, followed by synovitis, and progressing to OA as the clinical endpoint.² The demonstration that significant synovitis precedes structural changes in the progression of OA indicates that early intervention targeting pre-radiographic joint inflammation can delay or prevent chronic arthritic changes.³ Surgical and

non-surgical **synovectomy (synoviorthesis)** have been used for relieving synovitis in human, canine, and equine patients. This report describes **radiosynoviorthesis (RSO)** using a homogeneous **colloid** radiolabeled with **tin-117m**, a novel **radionuclide** that offers significant advantages over other radionuclides.



RADIOSYNOVIORTHESIS in clinical practice

The term “radiosynoviorthesis” (RSO) was introduced in Europe in the 1960s by Florian Delbarre to describe therapeutically active irradiation of the synovial lining.⁸ Rheumatologists administered a colloid embedded with a radionuclide (i.e., a **radiocolloid**) of yttrium-90 into the articular space. Using this process, the colloidal particles are phagocytized by macrophages in the synovial lining, after which they emit therapeutically active radiation within the synovial tissue until the radionuclide decays to its stable state. The number of inflammatory cells causing synovitis is reduced, and inflamed tissue is replaced with a fibrotic synovial membrane, with a corresponding alleviation of pain and improvement in function.^{9,10} An early study using beta-emitting samarium-153 in horses had mixed results, producing effective synovectomy but with transient lameness and swelling and exposure of some non-targeted periarticular tissue.¹¹

A key aspect of RSO is the choice of a radionuclide. Three radionuclides are widely used in clinical practice to treat synovitis: yttrium-90 (⁹⁰Y), rhenium-186 (¹⁸⁶Re), and erbium-169 (¹⁶⁹Er), all of which are artificially produced in a nuclear reactor.^{9,10,12}

In the case of **radiosynovectomy (RSV)** treatment, the radionuclide emits radiation that penetrates to the innermost layer of the synovium (synovial intima), where it produces energy of sufficient duration and intensity to achieve apoptosis and ablation of the inflamed cells. For this to occur, the radionuclide must have an adequate half-life ($t_{1/2}$), a selective tissue penetration range approximating the synovial thickness, and sufficient energy for therapeutic effect.

As ⁹⁰Y, ¹⁸⁶Re, and ¹⁶⁹Er decay, they emit radiation in the form of **beta particles** with a relatively wide tissue penetration range (Figure 1). While these radionuclides are therapeutically useful and have been evaluated in large clinical trials,⁹ their physical properties are not necessarily ideal for RSO. For example, ⁹⁰Y emits beta radiation that has a relatively wide range of soft tissue penetration, which risks irradiation of adjacent non-synovial tissue. ⁹⁰Y and ¹⁸⁶Re have short half-lives (2.7 and 3.7 days respectively), which creates storage and logistical limitations and may not consistently deliver sufficient radiation at the synovial target site.¹³

RADIOSYNOVIORTHESIS for intra-articular therapy

RSO has been successfully used in human medicine for more than 60 years in many countries, particularly in Europe where it was first described and where its use conforms to guidelines published by the European Association of Nuclear Medicine.^{8-10,12} RSO has been an accepted outpatient therapy for treatment of early stage chronic synovitis in rheumatoid arthritis, psoriatic arthritis, and OA patients for decades.^{10,14} RSO has important advantages versus surgical resection, the oldest ablative method. For example, as a minimally invasive procedure, RSO lowers bleeding risk in cases of hemophilic synovitis, where it is routinely used.^{13,15} As a localized treatment, RSO avoids problems associated with systemic therapies such as toxicity resulting from chronic use of non-steroidal anti-inflammatory drugs (NSAIDs) or immunosuppressive drugs. RSO also avoids tissue degradation that can occur from overuse of intra-articular corticosteroids. In human medicine, RSO has a favorable cost-benefit ratio, particularly when compared to surgery; a low rate of side effects; and application to virtually all articular joints, especially small, peripheral joints such as the finger joints.¹⁰ Current standards in human clinical practice generally take a conservative approach by recommending initial treatment with front-line therapies, including systemic NSAIDs,

glucocorticoids, and local joint therapies such as corticosteroid and hyaluronic acid injections prior to RSO.⁹ However, in patients that either respond poorly or have adverse side effects following these traditional therapies, RSO is a useful option that is now being considered in veterinary medicine.



TIN-117m: A NOVEL RADIONUCLIDE

Tin-117m is a unique radionuclide without the disadvantages of high-energy beta-emitting radionuclides. (Table 1 compares physical properties of tin-117m with other therapeutic radionuclides.)¹⁶ As such, tin-117m is particularly well suited for RSO in dogs. Instead of high-energy beta particles with a wide tissue penetration range (50–5,000 µm), tin-117m emits abundant **conversion electrons** (see Definitions), low-energy particles with a short penetration range of approximately

300 µm in tissue (Figure 1). Tin-117m has a $t_{1/2}$ of nearly 14 days, providing an ideal duration of effect spanning several half-lives to achieve therapeutic results and enable short-term stability during storage and handling. To illustrate, Figure 2 shows >99% dose retention in the joint of a dog one week following intra-articular injection with **homogeneous tin-117m colloid (HTC)**.¹⁷ No other radionuclide with the properties of tin-117m exists.¹⁸

Table 1. Comparison of radionuclides commonly used for radiosynoviorthesis¹⁶

Radionuclide	Half-life (days)	Avg. energy (maximum), keV	Avg. tissue penetration (maximum range), µm	Therapeutic emission	Diagnostic emission (keV)
Yttrium-90	2.7	934 (2,280)	4,100 (11,000)	beta	none
Rhenium-186	3.7	347 (1,070)	1,100 (4,400)	beta	gamma (137)
Erbium-169	9.4	100 (350)	140 (1,100)	beta	gamma (171)
Tin-117m	13.6	140 (152)	290 (300)	conversion electrons	gamma (159)

keV = kiloelectron volt.




Figure 1. The diagram compares the radiation dose range of conversion electrons emitted by tin-117m (300 μm , 0.3 mm) with beta radiation emitted by radionuclides such as yttrium-90 and erbium-169 (up to 11,000 μm , 11.0 mm). The ultra-narrow, discrete radiation range of tin-117m enables more precise dosimetry and avoidance of adverse effects on adjacent tissues such as what can occur with beta-emitting radionuclides.



Figure 2. Scintigraphy of an HTC-injected canine elbow shows high dose-retention of the homogeneous colloid with minimal uptake in the draining lymph node three days after administration. Retention at this time point was measured at >99% in synovial tissue, indicating a continuous therapeutic effect consistent with the 14-day half-life of tin-117m.

Image courtesy of Jimmy Lattimer, DVM, MS, DACVR

In addition to conversion electrons, tin-117m emits gamma radiation, a zero-mass quantum of light and electromagnetic radiation that results from nuclear decay of a radionuclide. Gamma radiation is nontherapeutic but readily detectable in tissue by imaging methods such as scintigraphy. By emitting gamma radiation at 159 kiloelectron volts (keV), tin-117m can be used diagnostically to detect the distribution and duration of its presence in tissue of treated patients. This application is similar to that for technetium-99m ($^{99\text{m}}\text{Tc}$), a widely used systemic radionuclide with gamma emissions of 140 keV that is used in diagnostic procedures, including evaluation of bone structure and function.

Due to its unique therapeutic and diagnostic (theranostic) properties as a conversion electron- and gamma-emitter with an optimal $t_{1/2}$, tin-117m has attracted interest as a radiopharmaceutical and also now as a medical device in the colloid form. Favorable results were reported in phase I and II clinical trials where tin-117m was used to treat metastatic bone pain in human patients.¹⁹⁻²¹ Investigators noted the value of the gamma emission component of tin-117m, which provides an objective basis for diagnostic monitoring, disease staging, dosage estimates, and assessing response to therapy.^{21,22}

A homogeneous colloid of tin-117m

Exubrion Therapeutics, Inc., has developed a patented preparation of tin-117m specifically for RSO and other potential applications in veterinary and human medicine. Tin-117m is manufactured using methods that produce yields sufficient to be scaled up for manufacturing therapeutic dosages in commercial quantities.¹⁸ The tin-117m radionuclide is combined with a homogeneous colloid.¹⁸ These colloidal objects clump together to form microparticles (1.5–20 µm) that are small enough to be phagocytized by synovial macrophages but large enough to avoid leakage outside the joint prior to phagocytosis. In situ retention of the HTC in laboratory animals has been measured out to five $t_{1/2}$ (i.e., 70 days), a duration sufficient for therapeutic effectiveness. The HTC has demonstrated safety and effectiveness following RSO of experimental OA in rats and safety in normal canine elbow joints (Figure 3).

Figure 3. Intra-articular injection of the radionuclide tin-117m into the caudolateral aspect of a canine elbow, positioned at 45-degree flexion, between the lateral condyle of the humerus and the triceps tendon. Following injection, the joint is put through a range of motion to disperse the radiocolloid throughout the synovial surface.

Image is courtesy of Caleb Hudson, DVM, DACVS, and Michelle Fabiani, DVM, DACVR, Gulf Coast Veterinary Specialists



CLINICALLY IMPORTANT FEATURES OF TIN-117M

Several features of tin-117m make it well suited for RSO and an improvement over other therapeutic radionuclides:

- **Localized administration:** Intra-articular dosing is suitable for outpatient use.
- **Non-beta emitter:** Avoids high-energy irradiation of non-synovial tissue, extra-articular diffusion, or systemic distribution.
- **Emits low-energy conversion electrons:** Minimizes potential for synovial scarring and eliminates collateral tissue damage.
- **Gamma radiation emitter:** Gamma energy of 159 keV is suitable for diagnostic imaging and is similar to the commonly used diagnostic radionuclide ^{99m}Tc (140 keV).
- **Half-life of 14 days:** Enables sufficient tissue retention for therapeutic effectiveness and a shelf life of 5 weeks.
- **Practical handling characteristics:** Ease of handling, hospital containment, and shipping using standard radiological safety and packaging practices.

RADIOSYNOVIORTHESIS in veterinary medicine

Radiotherapy has had various applications in companion animal medicine. For example, the beta emitter iodine-131 (^{131}I) has been used systemically to treat feline hyperthyroidism since the 1990s and is considered the treatment of choice for that condition.²³ Palliative and curative radiation therapy is now commonly used at veterinary oncology referral centers,²⁴ and radionuclides are also used for bone scanning in animals. Not surprisingly, successful RSO in human patients has created interest in using this method in companion animals, including horses, as a treatment for synovitis. Experimental RSO in horses has been attempted at university centers both in Europe and the U.S.^{11,25,26} Investigators in those studies used the beta-emitting radionuclides holmium-166 (^{166}Ho) or samarium-153 (^{153}Sm). However, high-energy emissions from either radionuclide resulted in some transient, periarticular soft tissue injury and minor extra-articular joint leakage.^{11,25,26} In a small Australian study, ^{90}Y was administered concurrently with methylprednisolone acetate to four horses with severe chronic synovitis and

hemarthrosis.²⁷ Median return to normal joint use was seven months, with two of the horses developing recurrent hemarthrosis.

In experimental studies with thulium-170 and ^{90}Y , healthy dogs were used as models for comparing results of canine and human RSO.²⁸⁻³⁰ Results indicated that RSO in dogs is feasible and generally well tolerated. However, the studies found that excessive dosages of beta-emitting radionuclides can reduce glycosaminoglycan synthesis in articular cartilage and result in extra-articular leakage of radiocolloid particles as late as nine months after intra-articular administration. Such outcomes reflect the importance in clinical applications of using well characterized radionuclides that emit radiation within well-defined parameters. Successful RSO in relatively small canine joints was noteworthy given the commonplace occurrence of canine elbow dysplasia and associated OA, a small-joint pathology that would be difficult to treat surgically.³¹

A safe and effective RSO radionuclide

Based on widespread clinical and experimental experience with beta-emitting radionuclides, a safe and effective radiocolloid suitable for RSO has the following characteristics:^{9,21}

- A limited, discrete emission penetration depth that corresponds to the thickness of the synovium (i.e., avoids irradiation that is too shallow for clinical effect or that extends beyond the synovial layer to affect non-target tissue).
- An intermediate radionuclide $t_{1/2}$ that is long enough to provide a reasonable shelf life and produce a therapeutic effect but short enough to avoid excessive exposure.
- A homogeneous colloid that binds the radionuclide so that it cannot escape beyond the joint.
- Suitable colloid particle size for synovial phagocytosis and in situ retention.

- Gamma emission for purposes of diagnostic imaging.
- High specific activity that allows scaling up for cost-effective production.
- A clinical profile that demonstrates a high degree of effectiveness and safety.

Because a non-beta-emitting, homogeneous tin-117m radiocolloid satisfies all of these criteria, it is considered to be uniquely suited for RSO treatment.¹⁸ Further evaluation in canine, feline, and equine models is expected to affirm its suitability for synovitis treatment combined with diagnostic confirmation of therapeutic response.

ACKNOWLEDGEMENTS

The authors acknowledge the contribution of Mark Dana of Scientific Communication Services, LLC, in the preparation of this report.

References

1. Benito MJ, Veale DJ, FitzGerald O, et al. Synovial tissue inflammation in early and late osteoarthritis. *Ann Rheum Dis*. 2005;64:1263-1267.
2. Sokolove J, Lepus CM. Role of inflammation in the pathogenesis of osteoarthritis: latest findings and interpretations. *Ther Adv Musculoskelet Dis*. 2013;5:77-94.
3. Scanzello CR, Umoh E, Pessler F, et al. Local cytokine profiles in knee osteoarthritis: elevated synovial fluid interleukin-15 differentiates early from end-stage disease. *Osteoarthritis Cartilage*. 2009;17:1040-1048.
4. Bondeson J, Wainwright SD, Lauder S, et al. The role of synovial macrophages and macrophage-produced cytokines in driving aggrecanases, matrix metalloproteinases, and other destructive and inflammatory responses in osteoarthritis. *Arthritis Res Ther*. 2006;8:R187.
5. de Lange-Brokaar BJ, Ioan-Facsinay A, van Osch GJ, et al. Synovial inflammation, immune cells and their cytokines in osteoarthritis: a review. *Osteoarthritis Cartilage*. 2012;12:1484-1499.
6. Poole AR. An introduction to the pathophysiology of osteoarthritis. *Front Biosci*. 1999;4:D662-D670.
7. McDougall JJ. Arthritis and pain: Neurogenic origin of joint pain. *Arthritis Res Ther*. 2006;8:220.
8. Delbarre F, Cayla J, Menkes C, et al. [Synoviorthesis with radioisotopes]. *Presse Med*. 1968;76:1045-1050.
9. Kampen WU, Voth M, Pinkert J, et al. Therapeutic status of radiosynoviorthesis of the knee with yttrium [⁹⁰Y] colloid in rheumatoid arthritis and related indications. *Rheumatology*. 2007;46:16-24.
10. Karavida N, Notopoulos A. Radiation synovectomy: an effective alternative treatment for inflamed small joints. *Hippokratia*. 2010;14:22-27.
11. Yarbrough TB, Lee MR, Hornof WJ, et al. Samarium 153-labeled hydroxyapatite microspheres for radiation synovectomy in the horse: a study of the biokinetics, dosimetry, clinical, and morphologic response in normal metacarpophalangeal and metatarsophalangeal joints. *Vet Surg*. 2000;29:191-199.
12. Klett R, Lange U, Haas H, et al. Radiosynoviorthesis of medium-sized joints with rhenium-186-sulfide colloid: a review of the literature. *Rheumatology*. 2007;46:1531-1537.
13. Silva M, Luck JV Jr, Llinas A. Chronic hemophilic synovitis: The role of radiosynovectomy. *Treatment Hemophilia*. 2004;33:1-10.
14. Modder G. Rheumatoid and related joint diseases. In: Radiosynoviorthesis. *Involvement of Nuclear Medicine in Rheumatology and Orthopaedics*. Meckenheim, Germany: Warlich Druck Verlagsges, MbH; 1995:13-23.
15. Rodriguez-Merchan EC, Wiedel JD. General principles and indications of synoviorthesis (medical synovectomy) in haemophilia. *Haemophilia*. 2001;7Suppl2:6-10.
16. Brenner W. Radionuclide joint therapy. In: Eary JF, Brenner W, eds. *Nuclear Medicine Therapy*. New York: Informa Healthcare; 2007:21-44.
17. Stevenson N, Lattimer J, Selting K, et al. Abstract S6-03: Homogeneous Sn-117m colloid—A novel radiosynovectomy agent. *World J Nucl Med*. 2015;14(Suppl 1):S15-S68.
18. Stevenson NR, St. George G, Simon J, et al. Methods of producing high specific activity Sn-117m with commercial cyclotrons. *J Radioanal Nucl Chem*. 2015;305:99-108.
19. Atkins HL, Mausner LF, Srivastava SC, et al. Tin-117m(4+)-DTPA for palliation of pain from osseous metastases: a pilot study. *J Nucl Med*. 1995;36:725-729.
20. Krishnamurthy GT, Swailem FM, Srivastava SC, et al. Tin-117m(4+)DTPA: pharmacokinetics and imaging characteristics in patients with metastatic bone pain. *J Nucl Med*. 1997;38:230-237.
21. Srivastava SC, Atkins HL, Krishnamurthy GT, et al. Treatment of metastatic bone pain with tin-117m Stannic diethylenetriaminepentaacetic acid: a phase I/II clinical study. *Clin Cancer Res*. 1998;4:61-68.
22. Srivastava SC. The role of electron-emitting radiopharmaceuticals in the palliative treatment of metastatic bone pain and for radiosynovectomy: applications of conversion electron emitter Tin-117m. *Brazilian Arch Biol Technol*. 2007;50:49-62.
23. Mooney CT. Radioactive iodine therapy in feline hyperthyroidism. *Vet Rec* 1990;127:555.
24. LaRue SM, Custis JT. Advances in veterinary radiation therapy: targeting tumors and improving patient comfort. *Vet Clin North Am Small Anim Pract*. 2014;44:909-923.
25. Mäkelä OT, Lammi MJ, Uusitalo H, et al. Effect of radiosynovectomy with holmium-166 ferric hydroxide macroaggregate on adult equine cartilage. *J Rheumatol*. 2004;31:321-328.
26. Mäkelä O, Sukura A, Penttilä P, et al. Radiation synovectomy with holmium-166 ferric hydroxide macroaggregate in equine metacarpophalangeal and metatarsophalangeal joints. *Vet Surg*. 2003;32:402-409.
27. Vallance SA, Lumsden JM, Begg Ap, et al. Idiopathic haemarthrosis in eight horses. *Aust Vet J*. 2012;90:214-220.
28. Hugenberg ST, Myers SL, Brandt KD. Suppression of glycosaminoglycan synthesis by articular cartilage, but not of hyaluronic acid synthesis by synovium, after exposure to radiation. *Arthritis Rheum*. 1989;32:4689-474.
29. Myers SL, Slowman SD, Brandt KD. Radiation synovectomy stimulates glycosaminoglycan synthesis by normal articular cartilage. *J Lab Clin Med*. 1989;114:27-35.
30. Polyak A, Das T, Chakraborty S, et al. Thulium-170-labeled microparticles for local radiotherapy: preliminary studies. *Cancer Biother Radiopharm*. 2014;29:330-338.
31. Kunst CM, Pease AP, Nelson NC, et al. Computed tomographic identification of dysplasia and progression of osteoarthritis in dog elbows previously assigned OFA grades 0 and 1. *Vet Radiol Ultrasound*. 2014;55:511-520.